ELECTRICAL MUSCLE STIMULATOR

cefar[®] myo²

USERS MANUAL

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USERS MANUAL

| | | | | | |
|--|----|------------------------------------|----|--|-----|
| PRECAUTIONARY MEASURES | 6 | 5.4 SPORT | 17 | 5.7 MASSAGE | 31 |
| 1. INTRODUCTION | 7 | Warm-up | 17 | Muscle Relaxation | 31 |
| 1.1. Electrical Muscle Stimulation (EMS) | 7 | Maximum Force | 17 | Lymphatic Massage | 31 |
| 2. THE CEFAR MYO 2 | 8 | Resistance Force | 17 | Massage Programs | 31 |
| 2.1 Functions | 9 | Explosive Force | 17 | 5.8 PAIN RELIEF | 32 |
| 2.2 How to use the stimulator | 9 | Endurance | 18 | TENS | 32 |
| 2.3 Changing sequence | 10 | Muscle Relaxation | 18 | TENS Muscle Relaxation | 32 |
| 2.4 Pausing the program | 10 | Capillarisation | 18 | Pain Relief - TENS Programs | 33 |
| 3. EFFECTIVE USE | 11 | Recovery | 18 | 5.9 REHABILITATION | 35 |
| 3.1 Electrode Placement | 11 | Active Recovery | 19 | Resuming Activity | 35 |
| Accessories | 11 | SPORT PROGRAMS Trunk | 19 | Rehabilitation Programs | 35 |
| 3.2 Body Position | 11 | SPORT PROGRAMS Upper Limbs | 21 | | |
| 4. PROGRAM STRUCTURE | 13 | SPORT PROGRAMS Lower Limbs | 23 | | |
| 4.1 Program Level | 13 | 5.5 FITNESS | 26 | 6. REPLACEMENT OF BATTERIES | 36 |
| 4.2 Program Sequences | 13 | Toning | 26 | 7. CARE INSTRUCTIONS | 36 |
| 4.3 Active Rest | 13 | Firming | 26 | 8. TROUBLESHOOTING | 36 |
| 4.4 Modulated Frequency | 13 | Sculpting | 26 | 9. TECHNICAL DATA | 37 |
| 4.5 Modulated Pulse Duration | 13 | Bodybuilding | 26 | 10. KEY TO THE STIMULATOR SYMBOLS | 37 |
| 5. STIMULATION TYPES AND PROGRAMS | 14 | FITNESS PROGRAMS Trunk | 27 | 11. GLOSSARY | 38 |
| 5.1 Stimulation Types | 14 | FITNESS PROGRAMS Upper Limbs | 28 | ELECTRODE PLACEMENT GUIDE | 115 |
| 5.2 Program Overview | 14 | FITNESS PROGRAMS Lower Limbs | 29 | | |
| 5.3 Program Guide | 16 | 5.6 AESTHETICS | 30 | | |
| Discover CEFAR MYO 2 | 16 | Lipolysis | 30 | | |
| | | Drainage | 30 | | |
| | | Aesthetic Programs | 30 | | |

PRECAUTIONARY MEASURES

Warning

- People with pacemakers, intra cardiac defibrillators or other active implants, pregnant women, and children must not use the CEFAR MYO 2.
- Do not stimulate on the front or sides of the neck, since a drop in blood pressure can occur.
- Do not use the CEFAR MYO 2 if you are connected to high-frequency surgical equipment, since this may cause burn injuries on the skin under the electrodes, as well as problems with the stimulator.
- Do not use the CEFAR MYO 2 in the direct vicinity of short-wave or microwave equipment (e.g. mobile phones), since this may affect the output power of the stimulator.
- Keep out of reach from children.

Note

- The stimulator should only be used with skin electrodes intended for nerve and muscle stimulation. See electrode package for instructions of use.
- The electrodes should only be placed on healthy skin. Avoid skin irritation by ensuring good contact between electrodes and skin. Hypersensitivity can occur in rare cases. If skin irritation should occur, treatment with the CEFAR MYO 2 should be temporarily discontinued. If condition persists, seek medical advice.

- Switch off stimulation before removing the electrodes from the skin. If an electrode comes off, shut off the stimulation before picking it up. Getting electrical stimulation through the fingers is unpleasant but not harmful.
- Exercise caution when using electrodes smaller than 14 cm², as the current density in some cases exceeds 2 mA rms/cm² and increases risk for burn injuries.
- Be careful with stimulating in the heart region.
- CEFAR Medical AB does not take any responsibility for any electrode placements other than recommended.
- Muscle soreness may occur after stimulation with CEFAR MYO 2, but usually disappears within a week.
- Only Cefar accessories should be used with the CEFAR MYO 2.
- Use the CEFAR MYO 2 only as stated in the operating instructions.
- Inspect the equipment prior to use.
- Never use more than one stimulator at the time.
- Never open the battery lid during stimulation.
- The stimulator meets the requirements of MDD 93/42/EEC, including the EMC requirements. This, however, is no guarantee that the equipment will not be affected by electromagnetic disturbances.
- If the stimulator is not used for some time (approximately 3 months), the non-rechargeable batteries should be removed from the stimulator.

1. INTRODUCTION

Thank you for choosing the CEFAR MYO 2. Your new stimulator has been designed with the latest technology in the field of electrical muscle stimulation. Using the CEFAR MYO 2 regularly will give your training an extra boost, regardless if you are a top athlete, or just exercise to keep fit.

In this manual you are shown how to operate your stimulator safely and how to use it in your training. Be sure to read it carefully before you start training to gain maximum benefits from the CEFAR MYO 2.

More information on our products and electrical muscle stimulation can be found on our web site: www.cefar.se

Good Luck!

1.1. Electrical Muscle Stimulation (EMS)


Electrical Muscle Stimulation (EMS) is successfully used in medical rehabilitation, and to support training in sports at all levels.

The main objective with EMS is to create contractions or vibrations in the muscles. Normally muscular activity is controlled by the central and peripheral nervous systems with electrical signals being sent to the muscles. EMS works in a similar way, but instead the muscle is stimulated with external electrical impulses acting through the skin via the nerve to the muscle. In both cases the muscle reacts and is contracted.

The main benefits of using EMS as a complement to your training

- During the electrical stimulation, all muscle fibres work simultaneously, making it a very efficient form of training compared to traditional muscle-building exercises.
- You can choose the exact muscle group you wish to train, without the fatigue and energy expenditure resulting from traditional muscle building exercises.
- EMS gives you a faster recovery since the blood circulation is increased and toxin elimination improved.
- EMS provides stimulation of the slow as well as the rapid muscle fibres. The rapid fibres are normally difficult to reach with ordinary training.
- Quick return to sport activities after injuries.

Arrow buttons ▲ ▼. Used for individually setting the amplitude for the different channels. **NOTE!** Increase the amplitude cautiously.

2.1 Functions

The CEFAR MYO 2 is a two-channel stimulator, i.e. one large group of muscles or two different parts of the body can be stimulated at the same time. The channels work simultaneously, i.e. both channels stimulate using the same program.

When using the CEFAR MYO 2 you run different programs designed for different types of stimulation. The stimulator has 67 programs organised in the following areas:

- Sport
- Fitness
- Aesthetics
- Massage
- Pain relief
- Rehabilitation

The CEFAR MYO 2 runs on two disposable 1.5 V AA or rechargeable 1.2 V AA batteries. A battery indicator shows when the batteries need to be replaced. After using the stimulator it automatically shuts off after 5 minutes of inactivity to save the batteries.

The stimulator has a belt clip that enables you to have your hands free during stimulation.

2.2 How to use the stimulator

IMPORTANT! The amplitude can be increased and decreased at any time during stimulation using the ▲ and ▼ buttons. If you ever feel discomfort or pain immediately decrease the amplitude. You can also use the On/Off button ● to terminate the stimulation at all times.

1. Connect a pair of electrodes to a cable, and the cable to a channel on the stimulator.

2. Refer to the electrode placement guide in the back of this manual to see where the electrodes should be placed. Wash and dry the area of the skin where you will be placing the electrodes. See section 3 for more information on electrode placement.
3. Attach the electrodes to your body.
4. Switch the stimulator on by pressing the  button. A beep will indicate that it is switched on. The stimulator automatically chooses the program used the previous time. If you wish to change the program, press the **P+** or **P-** buttons to browse through the different choices. If you are using your CEFAR MYO 2 for the first time we suggest you use Program 0, which is a test program. See section 5.3 for more information on this program.

5. Start the stimulation by setting the desired amplitude with ▲ and ▼ for each channel. The amplitude is set in steps of 0.5 mA which allows you to adjust it exactly after your individual preferences. Each channel is set separately. Note: With intermittent stimulation an amplitude for the rest time needs to be set if the rest is active, see section 4.3 for more information on Active Rest.
6. When the program is finished, three beeps will sound and the timer symbol in the display will blink. The amplitude will decrease to 0.0 mA automatically.
7. You can always stop stimulation by using the  button at all times.

A program is made up of 1, 2 or 3 sequences, all with different purposes. It could be a warm-up sequence followed by a training sequence and a recovery sequence. A beep is heard every time a program changes sequences. When one sequence ends, another starts immediately. You can change sequences whenever you want by pressing the **P+** button or **P-** button for 5 seconds. When the program is changing sequence the amplitude drops to half the setting to prevent discomfort during the transition. This means that you have to raise the amplitude at the start of each new sequence.

You can pause the program by setting the amplitude to 0.0 mA. Start stimulation again by increasing the amplitude. The program will resume at the point at which stimulation was interrupted.

PROGRAM

| BODY PART | Maximum Force, Resistance Force, Explosive Force, Endurance, Toning, Firming, Sculpting, Bodybuilding, Resuming Activity | Warm-up, Muscle Relaxation, Massage programs, Drainage, Pain Relief programs, Lipolysis, Capillarisation, Recovery, Active Recovery |
|---|--|---|
|  <p>QUADRICEPS</p>  <p>LOWER FRONT LEG</p>  <p>CALF</p>  <p>KNEE *</p>  <p>ANKLE *</p> |  |   |
|  <p>ABDOMINAL</p>  <p>HIPS</p> |  |   |
|  <p>LUMBAR/BACK</p>  <p>LATISSIMUS DORSI</p> |   |   |
|  <p>HAMSTRINGS</p>  <p>GLUTEAL/BUTTOCKS</p> |  |   |
|  <p>TRICEPS</p>  <p>SHOULDER</p>  <p>NECK *</p> |  |   |
|  <p>BICEPS</p>  <p>FOREARM</p>  <p>WRIST *</p>  <p>ELBOW *</p> |  |   |
|  <p>PECTORAL/CHEST</p> |  |   |

* PAIN RELIEF

4. PROGRAM STRUCTURE

When using the CEFAR MYO 2 you run different programs with different types of stimulation. For example, the programs for endurance training are totally different from programs for after-training recovery. Furthermore, some of the programs have different levels. The programs are also adapted for different body parts.

4.1 Program Level

Some of the programs have three different levels: Level 1, 2 and 3, e.g. Maximum Force 1, Maximum Force 2 and Maximum Force 3. Use these levels to gradually increase your training load. In the programs with different levels, each level is recommended to be used 6-8 weeks before going on to the next level.

4.2 Program Sequences

Every program consists of 1-3 sequences, each with its own purpose. The programs with 3 sequences always start with a warm-up sequence and end with a recovery sequence to ensure a balanced workout.

| | | |
|------------|----------------------------------|----------------------------|
| Sequence 1 | Warm-up | Moderate muscle vibrations |
| Sequence 2 | ----- See Stimulation Type ----- | |
| Sequence 3 | Recovery | Visible muscle vibrations |

When the program is changing sequences, the amplitude drops to half its value for a more comfortable transition. This means you have to raise the amplitude at the start of each new sequence.

On the display of the stimulator you can see how many sequences the program consists of (number of squares), the current sequence (the filled square), and the duration of the sequence. You may change the stimulation time for the current sequence, or switch to the next or previous sequence at any time during stimulation.

4.3 Active Rest

In the training programs where the intensity is high enough for muscle contractions, there is stimulation between the contractions, so-called Active Rest. The Active Rest stimulation should cause muscle vibrations, but no contractions. The aim is to keep the muscle prepared for the next contraction, and to decrease muscle soreness afterwards.

You have to raise the amplitude for the Active Rest $\cdots \wedge \cdots$ as well as for the contractions $\text{---} \nearrow \text{---} \searrow \text{---}$.

4.4 Modulated Frequency

Modulated Frequency means that the frequency varies over time, and can in some cases make the stimulation more effective and comfortable. This type of stimulation is used in programs such as Sculpting and Maximum Force.

4.5 Modulated Pulse Duration

This type of stimulation is used to make the stimulation more comfortable. The width of each pulse varies over a certain time. The feeling of this stimulation type is "wave-like" and used e.g. in the Muscle Relaxation programs.

The programs in the CEFAR MYO 2 are categorised in different application areas: Sport, Fitness, Aesthetics, Massage, Pain Relief/TENS and Rehabilitation. Every area has different types of stimulation.

5.1 STIMULATION TYPES

| | |
|-------------------|-------------------------------|
| SPORT | AESTHETICS |
| Warm-up | Lipolysis |
| Maximum Force | Drainage |
| Resistance Force | |
| Explosive Force | MASSAGE |
| Endurance | Muscle Relaxation |
| Muscle Relaxation | Muscle Relaxation Face |
| Capillarisation | Lymphatic Massage |
| Recovery | |
| Active Recovery | PAIN RELIEF |
| | TENS |
| | TENS Muscle Relaxation |
| FITNESS | |
| Toning | |
| Firming | REHABILITATION |
| Sculpting | Resuming Activity |
| Bodybuilding | |

TRUNK

UPPER LIMBS

LOWER LIMBS

5.2 PROGRAM OVERVIEW

| NO. 0 | BODYPART | TITLE Test Program |
|------------------------|----------|-----------------------|
| SPORT - TRUNK | | |
| 1 | Sport | Warm-up |
| 2 | Sport | Maximum Force |
| 3 | Sport | Resistance Force |
| 4 | Sport | Explosive Force |
| 5 | Sport | Endurance |
| 6 | Sport | Muscle Relaxation |
| 7 | Sport | Recovery |
| 8 | Sport | Active Recovery |
| SPORT - UPPER LIMBS | | |
| 9 | Sport | Warm-up |
| 10 | Sport | Maximum Force 1 |
| 11 | Sport | Maximum Force 2 |
| 12 | Sport | Maximum Force 3 |
| 13 | Sport | Resistance Force 1 |
| 14 | Sport | Resistance Force 2 |

| | | |
|----|-------|-------------------|
| 15 | Sport | Explosive Force 1 |
| 16 | Sport | Explosive Force 2 |
| 17 | Sport | Endurance |
| 18 | Sport | Muscle Relaxation |
| 19 | Sport | Capillarisation |
| 20 | Sport | Recovery |
| 21 | Sport | Active Recovery |

SPORT

- LOWER LIMBS

| SPORT - LOWER LIMBS | | |
|---------------------|-------|--------------------|
| 22 | Sport | Warm-up |
| 23 | Sport | Maximum Force 1 |
| 24 | Sport | Maximum Force 2 |
| 25 | Sport | Maximum Force 3 |
| 26 | Sport | Resistance Force 1 |
| 27 | Sport | Resistance Force 2 |
| 28 | Sport | Resistance Force 3 |
| 29 | Sport | Explosive Force 1 |
| 30 | Sport | Explosive Force 2 |
| 31 | Sport | Explosive Force 3 |
| 32 | Sport | Endurance 1 |
| 33 | Sport | Endurance 2 |
| 34 | Sport | Muscle Relaxation |
| 35 | Sport | Capillarisation |
| 36 | Sport | Recovery |
| 37 | Sport | Active Recovery |

FITNESS

- TRUNK

| FITNESS | | |
|---------|---------|--------------|
| - TRUNK | | |
| 38 | Fitness | Toning |
| 39 | Fitness | Firming |
| 40 | Fitness | Sculpting |
| 41 | Fitness | Bodybuilding |

FITNESS

- UPPER LIMBS

| FITNESS | | |
|---------------|---------|--------------|
| - UPPER LIMBS | | |
| 42 | Fitness | Toning |
| 43 | Fitness | Firming |
| 44 | Fitness | Sculpting |
| 45 | Fitness | Bodybuilding |

FITNESS

- LOWER LIMBS

| FITNESS | | |
|---------------|---------|--------------|
| - LOWER LIMBS | | |
| 46 | Fitness | Toning |
| 47 | Fitness | Firming |
| 48 | Fitness | Sculpting |
| 49 | Fitness | Bodybuilding |

AESTHETICS

| AESTHETICS | | |
|------------|------------|-----------|
| 50 | Aesthetics | Lipolysis |
| 51 | Aesthetics | Drainage |

MASSAGE

| MASSAGE | | |
|---------|---------|------------------------|
| 52 | Massage | Muscle Relaxation |
| 53 | Massage | Muscle Relaxation Face |
| 54 | Massage | Lymphatic Massage |

PAIN RELIEF

| PAIN RELIEF | | |
|-------------|------|--|
| 55 | TENS | In general |
| 56 | TENS | Neck |
| 57 | TENS | Shoulder |
| 58 | TENS | Elbow |
| 59 | TENS | Lumbar |
| 60 | TENS | Knee |
| 61 | TENS | Ankle |
| 62 | TENS | TENS Muscle Relaxation, Upper limbs |
| 63 | TENS | TENS Muscle Relaxation, Lower limbs |
| 64 | TENS | TENS Muscle Relaxation, Trunk |

REHABILITATION

| REHABILITATION | | |
|----------------|----------------|--------------------------------|
| 65 | Rehabilitation | Resuming Activity, Trunk |
| 66 | Rehabilitation | Resuming Activity, Upper limbs |
| 67 | Rehabilitation | Resuming Activity, Lower limbs |

The Program Guide helps you to choose a program suitable for your individual needs. Each section specifies the following:

PROGRAM NUMBER

| | |
|----------------------------|--|
| APPLICATION AREA | Stimulation type |
| BODY PART/MUSCLE | Which part of the body to stimulate |
| SEQUENCE | Total duration of each sequence |
| TOTAL DURATION | Total duration of the program |
| INTENSITY | The strength of the stimulation |
| ELECTRODE PLACEMENT | A suggestion of electrode placement (see Electrode Placement Guide , page 115) |

Discover CEFAR MYO 2

If you are using the CEFAR MYO 2 for the first time, we suggest that you try PROGRAM 0 to familiarise yourself with the stimulator.

| | |
|----------------------------|------------------------------|
| PROGRAM 0. | TEST PROGRAM |
| Sequence | 2 min / 2 min / 2 min |
| Total duration | 6 minutes |
| Electrode placement | 11a |

Sequence 1 (duration 2 min)

Progressively increase the intensity level until you reach a visible vibration. You will feel that the vibrating sensation changes into a light contraction, then decreases to a vibration. This is continuous stimulation with modulated frequency. You can find this form of stimulation in e.g. the warm-up sequences, the Active Recovery programs and the Drainage programs.

Sequence 2 (duration 2 min)

You must increase the level of intensity again. You will feel muscle contractions of medium intensity, alternating with rest periods. The rest time is active, i.e. stimulation is active also during the rest time. The higher the intensity, the larger the number of stimulated muscle fibres. This is the most classical form of muscle stimulation. You will find this form of stimulation in the training programs e.g. in the Maximum Force programs.

Sequence 3 (duration 2 min)

You must increase the level of intensity again. This is continuous stimulation with modulated pulse duration. The feeling of this stimulation is "wave like" and used e. g in the Muscle Relaxation programs.

5.4 SPORT

Warm-up

Increases the muscle's temperature to make it more prepared for exercise. It also increases the blood supply and the metabolism in the muscle.

-  Recommended before physical training or competition.

Maximum Force

Improves the maximum-force capacity and increases muscle mass. Allows training with a maximum load with less risk of injuries compared to classic heavy strength training.

-  Recommended for all sports where you need maximum strength – 1-3 times per week.

Resistance Force

Improves the capacity to develop a high level of force over a prolonged period of time and increases the ability to withstand toxin accumulation and therefore delays the fatigue of the muscle. To achieve a very high intensity of work the objective is to induce a large number of long strong contractions alternating with short Active Rest.

- 1** This is recommended for sports needing intense muscle work over a prolonged period of time – 1-3 times per week.

Explosive Force

Explosive Force increases the capacity for instantaneous maximum force. The aim is to utilise the largest amount of muscle fibres in the shortest possible time. This type of program is used as a more comfortable complement to general explosive force training, which usually is very hard and demanding. It is recommended to start with a Maximum Force program to prepare the muscle for the Explosive Force.

-  The Explosive Force programs are used for sports where strength at a high speed in the shortest possible time is important e.g. sprint, ball sports, alpine skiing – 1-3 times per week.

Active Recovery

Increases the ability to regain muscular capacity after training or competition and to resume training more easily. The stimulation varies from muscle vibrations to light contractions. This increases circulation and decreases muscle stiffness. It also gives pain relief due to the creation of endorphins. This program is valuable in all sports as it improves the ability to return quickly to maximum muscular capacity.

 This program is recommended to be used 6-24 hours after training or competition.

SPORT PROGRAMS

Trunk

PROGRAM 1.

SPORT

Sequence

Total duration

Intensity

Electrode placement

WARM UP

5 min

5 minutes

Moderate muscle vibrations.

1a, 2a, 2d, 4a, 5a, 6a

PROGRAM 2.

SPORT

Sequence

Total duration

Intensity

Electrode placement

MAXIMUM FORCE

5 min / 16 min / 10 min

31 minutes

Maximum muscle contraction without pain.

The Active Rest should give muscle vibrations.

2a, 2d, 2e, 3a, 4a, 5a, 6a

PROGRAM 3.

SPORT

Sequence

Total duration

Intensity

Electrode placement

RESISTANCE FORCE

4 min / 14 min / 15 min

33 minutes

Medium to strong muscle contractions.

The Active Rest should give muscle vibrations.

2a, 2d, 2e, 3a, 4a, 5a, 6a

SPORT

Total duration

Intensity

Electrode placement

EXPLOSIVE FORCE

35 minutes

Maximum muscle contractions without pain.

The Active Rest should give muscle vibrations.

2a, 2e, 3a, 4a, 6a

SPORT

Total duration

Intensity

Electrode placement

ENDURANCE

60 minutes

Medium, visible muscle contractions.

The Active Rest should give muscle vibrations.

2a. 2d. 2e. 3a. 4a. 5a. 6a

SPORT

Total duration

Intensity

Electrode placement

MUSCLE RELAXATION

20 minutes

Visible muscle vibrations to light contractions.

1a, 2a, 2d, 2e, 3a, 4a, 5a, 6a

SPORT

Total duration

Intensity

Electrode placement

RECOVERY

10 minutes

Visible muscle vibrations.

1a, 2a, 2d, 2e, 3a, 4a, 5a, 6a

SPORT

Total duration

Intensity

Electrode placement

ACTIVE RECOVERY

20 minutes

Muscle vibrations and light contractions.

1a, 2a, 2d, 2e, 3a, 4a, 5a, 6a

Electrode placement

5 min / 15 min / 15 min
35 minutes

**Maximum muscle contractions without pain.
The Active Rest should give muscle vibrations.
7a, 7c, 7d, 9a, 9d, 11a, 11b**

Electrode placement

5 min / 40 min / 10 min

55 minutes
Medium, visible muscle contractions.
The Active Rest should give muscle vibrations.
7a, 7c, 7d, 9a, 9d, 11a, 11b

Sequence
Total duration
Intensity
Electrode plac

20 min

20 minutes
Visible muscle vibrations to light contractions.
7a, 7b, 9a, 9c, 9d, 11a, 11b, 11c

Sequence
Total duration
Intensity
Electrode placement

25 min

25 minutes
Visible muscle vibrations.
7a. 7b. 9a. 9c. 9d. 11a. 11b. 11c

Sequence
Total duration
Intensity
Electrode placement

10 min

10 minutes
Visible muscle vibrations.
7a, 7b, 9a, 9c, 9d, 11a, 11b, 11c

Sequence
Total duration
Intensity
Electrode placement

20 min

20 minutes
Muscle vibrations and light contractions.
7a, 7b, 9a, 9c, 9d, 11a, 11b, 11c

Electrode placement

5 min / 15 min / 15 min
35 minutes

**Maximum muscle contractions without pain.
The Active Rest should give muscle vibrations.
13a, 13d, 15a, 16a**

Electrode placement

5 min / 15 min / 15 min
35 minutes

**Maximum muscle contractions without pain.
The Active Rest should give muscle vibrations.
13a. 13d. 15a. 16a**

Electrode placement

5 min / 15 min / 15 min
35 minutes

**Maximum muscle contractions without pain.
The Active Rest should give muscle vibrations.
13a, 13d, 15a, 16a**

Electrode placement

5 min / 40 min / 10 min
55 minutes

Medium, visible muscle contractions.
The Active Rest should give muscle vibrations.
13a, 13d, 15a, 16a, 17a, 19a

Electrode placement

5 min / 40 min / 15 min
60 minutes

Medium, visible muscle contractions.
The Active Rest should give muscle vibrations.
13a. 13d. 15a. 16a. 17a. 19a

Electrode placement

20 min

20 minutes
Visible muscle vibrations to light contractions.
13a, 13c, 15a, 15b, 16a, 16b, 17a, 19a

PROGRAM 35.

SPORT

Sequence

Total duration

Intensity

Electrode placement

CAPILLARISATION

25 min

25 minutes

Visible muscle vibrations.

13a, 13c, 15a, 15b, 16a, 16b, 17a, 19a

PROGRAM 36.

SPORT

Sequence

Total duration

Intensity

Electrode placement

RECOVERY

10 min

10 minutes

Visible muscle vibrations.

13a, 13c, 15a, 15b, 16a, 16b, 17a, 19a

PROGRAM 37.

SPORT

Sequence

Total duration

Intensity

Electrode placement

ACTIVE RECOVERY

20 min

20 minutes

Muscle vibrations and light contractions.

13a, 13c, 15a, 15b, 16a, 16b, 17a, 19a

The toning programs offer genuine maintenance work for muscular tone. This type of training session is of moderate intensity, making it suitable for everyone, without causing muscle soreness.

- These programs represent training exercise of moderate intensity, aimed at improving muscular density without causing muscle soreness.

- i** Recommended for initiating the build up of weak muscles that have lost some of their bulk. They can be used every day on large muscle groups such as the thighs, abdominal and dorsal muscles.

The sculpting programs produce intense work comparable with that afforded by a regular strength training session. It is normal for some muscle soreness to be felt following sessions of this nature.

- The frequency of the impulses pushed to their maximum produces highly intense muscular work. This form of training is aimed at achieving the hypertrophy of the fast muscle fibres. It brings about complete muscular exhaustion, followed by intense muscle soreness.

-  It is aimed at those who are used to muscle training and wish to improve their muscular definition. These programs are recommended to be used 1-3 times per week.

FITNESS PROGRAMS

Upper Limbs

FITNESS

TONING

5 min / 15 min / 5 min

25 minutes

Medium muscle contractions.

7a, 7b, 7c, 9b, 9d, 11a, 11b

FITNESS

FIRMING

5 min / 20 min / 5 min

30 minutes

Maximum muscle contractions without pain.

The Active Rest should give muscle vibrations.

7a. 7b. 7c. 9b. 9d. 11a. 11b

FITNESS

SCULPTING

5 min / 20 min / 5 min

30 minutes

Maximum muscle contractions without pain.

The Active Rest should give muscle vibrations.

7a. 7b. 7c. 9b. 9d. 11a. 11b

FITNESS

BODYBUILDING

5 min / 15 min / 5 min

25 minutes

Maximum muscle contractions without pain.

The Active Rest should give muscle vibrations.

7a, 7b, 7c, 9b, 9d, 11a, 11b

FITNESS PROGRAMS

Lower Limbs

FITNESS

TONING

5 min / 15 min / 5 min

25 minutes

Medium muscle contractions.

13a, 15a, 16b, 17a, 19a

FITNESS

FIRMING

5 min / 20 min / 5 min

30 minutes

Maximum muscle contractions without pain.

The Active Rest should give muscle vibrations.

13a. 15a. 16b. 17a. 19a

FITNESS

SCULPTING

5 min / 20 min / 5 min

30 minutes

Maximum muscle contractions without pain.

The Active Rest should give muscle vibrations.

13a, 15a, 16b, 17a, 19a

FITNESS

BODYBUILDING

5 min / 15 min / 5 min

25 minutes

Maximum muscle contractions without pain.

The Active Rest should give muscle vibrations.

13a, 15a, 16b, 17a, 19a

The low frequency stimulation in this program has a positive effect on the metabolism, promoting the elimination of toxins and unwanted fluids. It can be used for reshaping and rejuvenating areas by modifying the metabolism of lipocytes. For best result, follow up with aerobic training e.g. running or cycling for a minimum of 20 minutes.

- ## Drainage

The co-ordinated stimulation achieved by the two channels of the stimulator produces a muscular pumping effect. The pumping action mainly activates the venous but also the lymphatic systems mechanically.

-  The objective is to improve the circulation and to promote the elimination of toxins.
Can be used every day.

Aesthetic Programs

Electrode placement

5b, 6a, 13a, 15a, 15b, 17a, 19a

Electrode placement

5b, 6a, 13a, 15a, 15b, 17a, 19a

5.7 MESSAGE

Muscle Relaxation

Enables a very effective decontraction of the muscle due to the comfortable vibrations that increase the circulation and help the muscles to relax.

 These programs can be used on any occasion when muscle relaxation is needed.

Lymphatic Massage

The coordinated stimulation achieved by the two channels of the stimulator produces a muscular pumping effect. The pumping action mainly activates the lymphatic but also the venous systems mechanically.

1 Its objective is to improve the circulation in the lymphatic system and to promote the elimination of toxins. These programs can be used on any occasion when Lymphatic Massage is needed.

PROGRAM 52.

MASSAGE

Sequence

Total duration

Intensity

Electrode placement

MUSCLE RELAXATION

20 min

20 minutes

Visible muscle vibrations.

**1a, 1b, 2a, 2d, 4b, 6a, 7a, 9c, 11c,
13a, 15a, 15b, 16b, 17a, 19a**

PROGRAM 53.

MASSAGE

Sequence

Total duration

Intensity

Electrode placement

MUSCLE RELAXATION FACE

15 min

15 minutes

Light muscle vibrations.

20a

PROGRAM 54.

MASSAGE

Sequence

Total duration

Intensity

Electrode placement

LYMPHATIC MASSAGE

20 min

20 minutes

Visible muscle vibrations to light contractions.

6a, 9c, 11c, 13a, 15a, 15b, 16b, 17a, 19a

TENS

TENS (Transcutaneous Electrical Nerve Stimulation) utilises the nervous system's own pain relief mechanisms to ease both acute and chronic pain. When stimulating, impulses are sent through the nerves to block the pain impulses. The pain relief is usually most efficient during stimulation, but the effect can also last after the program has finished. TENS also increases the circulation.

-  The TENS programs do not always treat the cause of pain. If the pain persists, seek medical advice.
-  Place the electrodes on the painful area. Every stimulation session should last at least 30 minutes and the program can be used on any occasion when pain relief is needed.

TENS Muscle Relaxation

Offers effective relaxation of the muscle, as well as pain relief. The TENS utilises the nervous system's own pain relief mechanisms to ease both acute and chronic pain. When stimulating, impulses are sent through the nerves to block the pain impulses. The muscle relaxation increases the circulation and helps the muscles to relax.

-  The TENS programs do not always treat the cause of pain. If the pain persists, seek medical advice.
-  As tensed muscles sometimes are painful this program is combined TENS and muscle relaxation. The first phase is TENS and the second is muscle relaxation. The program can be used on any occasion when both pain relief and muscle relaxation are needed.

Sequence
Total duration
Intensity

Electrode placement

25 min / 15 min
40 minutes

- **TENS sequence one:** Pins-and-needles sensation without pain. Should not cause contractions.
- **Muscle relaxation sequence two:** Visible muscle vibrations.

Sequence
Total duration
Intensity

Electrode placement

25 min / 15 min
40 minutes

- **TENS sequence one: Pins-and-needles sensation without pain. Should not cause contractions.**
- **Muscle relaxation sequence two: Visible muscle vibrations.**

Sequence
Total duration
Intensity

Electrode placement

25 min / 15 min
40 minutes

- **TENS sequence one:** Pins-and-needles sensation without pain. Should not cause contractions.
- **Muscle relaxation sequence two:** Visible muscle vibrations.

5.9 REHABILITATION

Resuming Activity

Resuming physical activity after a period of interruptions. After injuries or other interruption in physical activity the muscles often lose strength and muscle mass. It is important to restore normal muscle capacity as soon as possible. The Resuming Activity program enables you to train the specific muscles that are affected.

1 Can be used on any occasions when resuming activity is needed.

Rehabilitation Programs

PROGRAM 65. REHABILITATION

RESUMING ACTIVITY - TRUNK

| | |
|---------------------|-------------------------------|
| Sequence | 5 min / 30 min / 10 min |
| Total duration | 45 minutes |
| Intensity | Moderate muscle contractions. |
| Electrode placement | 2a, 2d, 2e, 3a, 4a, 5a, 6a |

PROGRAM 66. REHABILITATION

RESUMING ACTIVITY - UPPER LIMBS

| | |
|---------------------|-------------------------------|
| Sequence | 5 min / 30 min / 10 min |
| Total duration | 45 minutes |
| Intensity | Moderate muscle contractions. |
| Electrode placement | 7a, 7c, 7d, 9a, 9d, 11a, 11b |

PROGRAM 67. REHABILITATION

RESUMING ACTIVITY - LOWER LIMBS

| | |
|---------------------|-------------------------------|
| Sequence | 5 min / 30 min / 10 min |
| Total duration | 45 minutes |
| Intensity | Moderate muscle contractions. |
| Electrode placement | 13a, 13d, 15a, 16b, 17a, 19a |

Switch off the stimulator. Remove the battery cover by sliding it downwards and then lifting it upwards, and take out the batteries. Make sure the new batteries are placed correctly with respect to the polarity (note + and - symbols). See fig. page 113. Dispose the used batteries at a recycling station.

7. CARE INSTRUCTIONS

Taking care of and cleaning the CEFAR MYO 2 equipment is simple, provided the following instructions are observed:

- Always keep the stimulator and accessories in the original case when not in use.
- Never expose the stimulator to water. Wipe it off with a damp cloth if necessary.
- The leads are best preserved if left attached to the stimulator between sessions.
- Do not jerk cords or connections.

8. TROUBLESHOOTING

Check that all settings are correct and make sure that the electrodes are correctly placed (see the Electrode Placement Guide p. 115).

- The skin is irritated. For advice on skin care, see section Precautionary Measures.
- Electrodes are too old and need to be replaced.

Current pulses are weak, or missing

- Check if the battery needs replacing, see section 6, Replacement of batteries.
- Make sure that the electrodes do not touch each other.
- If the circuit is not closed, e.g. due to a broken cable, the amplitude for this channel is automatically turned off and 0.0 mA starts blinking. A cable breakage can be checked by pressing the cable's pins against one another while increasing the amplitude for this channel to about 15 mA. If the amplitude now turns off and 0.0 mA blinks, the cable needs to be replaced. NOTE! When checking for cable breakage the amplitude must not be increased to more than 20 mA, since this could cause damage to the stimulator.
- Excessively high resistance is indicated by the amplitude dropping to 0 mA and 0.0 mA blinking.
- The electrodes are worn out and need replacement.


Service of the CEFAR MYO 2 should only be performed by CEFAR Medical AB or Cefar distributors, and should be carried out once a year.

9. TECHNICAL DATA

CEFAR MYO 2 is an electrical stimulator, with two channels, intended for muscular training and pain relief. The stimulator runs on either two rechargeable 1.2 V AA batteries or two non-rechargeable 1.5 V AA batteries. It delivers a constant current of 100 mA through a load of 1000 ohm.

| | |
|-------------------------------|---|
| NUMBER OF CHANNELS | 2 (simultaneous) |
| CONSTANT CURRENT | Up to a resistance of 1000 ohm |
| STIMULATION FORMS | Continuous stimulation: constant pulse, modulated pulse duration, modulated frequency. Intermittent stimulation: constant pulse, modulated frequency. |
| AMPLITUDE | 0-100 mA, in steps of 0.5 mA. |
| WAVEFORM | Symmetrical biphasic square pulse |
| MAXIMUM FREQUENCY | 120 Hz |
| MAXIMUM PULSE DURATION | 300 μ s |
| POWER SOURCE | 2 x 1.5 V AA non-rechargeable or 2 x 1.2 V AA rechargeable batteries |

**CURRENT CONSUMPTION
FOR ONE CHANNEL, 300 μ S,
30 HZ, 30MA**

75 mA

ENVIRONMENT FOR STORAGE AND SHIPPING

Temperature 10° C - 40° C
Air humidity 30% - 75%
Air pressure 700 hPa-1060 hPa

EXTERNAL DIMENSIONS

65 x 120 x 30 mm

WEIGHT

approximately 160 g

I R.M.S. MAX/CHANNEL

27 mA

10. KEY TO THE STIMULATOR SYMBOLS

Read the operating instructions prior to use.


Patient part type – Body Floating.

Complies with MDD 93/42/EEC requirements.
Notified body SEMKO AB.

Warranty: The stimulator is covered by a two-year warranty against manufacturing faults. Subject to normal conditions of use.

| | |
|--|--|
| METABOLISM | Chemical processes that are continually taking place in the human body and which are essential to life e.g. the building of tissues and the elimination of waste matters. |
| MODULATED FREQUENCY | The frequency varies over a certain time. |
| MODULATED PULSE | |
| DURATION | The duration of each pulse varies over a certain time. |
| MUSCLE FIBRES, — SLOW, FAST | Component in the muscle. Muscles are made up of slow muscle fibres and fast muscle fibres. The slow fibres are mainly used in endurance sports, whereas the fast fibres are mainly used in explosive sports. |
| PULSE DURATION | The phase duration of each individual pulse. |
| SEQUENCE | Part of a program. |
| TENS | Transcutaneous Electrical Nerve Stimulation. |
| TOXINS | Poisonous substances in the body. |
| VENOUS SYSTEM | System of veins which bring blood back to the heart from the tissues. |
| Y-CABLE  | An electrode cable making it possible to double the amount of electrodes. |

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Photo © Studio Ekvall
MDD 93/42/EEC. Organisme notifié SEMKO AB.
Helsingborg AB
CEFAR Medical AB | Ideon Science Park